


# Special extruder BGEMH 40



(<http://rozrobkinauchni.vsocorp.org/>). E-mail: kbgarmasch@ukr.net

mob. tel. +38 095 466 23 89; +38 096 074 06 49



# SPECIAL EXTRUDER BGEMH 40

- ✦ Special molding head for the production of artificial protein shell of caliber 40. This device is similar to the scheme of the German engineering office of Dr. Dorstewitz.
- ✦ The peculiarity of the device is the ability to work on the classic scheme used at the Ukrainian enterprise LLC «Belkozin», as well as the axial supply of protein collagen. Details for interested parties.
- ✦ You can order a special extruder for the formation of an artificial protein shell by these phones or e-mail.
- ✦ E-mail: [kbgarmasch@ukr.net](mailto:kbgarmasch@ukr.net)  
mob. tel.: +38 095 466 23 89;
- ✦ +38 096 074 06 49





# SPECIAL EXTRUDER BGEMH 40

Caliber of the formed artificial protein shell, mm	40
The number of versions of the extruder	9
Extruder type	Rotary-toothed with a fixed shaft
The size of the gap in the area of weaving of the protein shell, mm	0,5
The ratio of the diameters of the rotor and the fixed shaft in the area of weaving the shell, mm/mm	$\varnothing 40H10/\varnothing 39h10$
The size of the slit of the calibration die, mm	$0,525 \pm 0,01$ $0,55 \pm 0,01$ $0,57 \pm 0,01$
Power of the driving electric motor, kW	1,5
Mains voltage, V	230/400
Power supply voltage frequency, Hz	50
Unit weight, kg	97

